

EN:14683:2019 Standard

For avoiding spreading and projections of the wearer, projection of secretions of saliva which may contain transmissible infectious agents:

- via droplets (transmission via droplets of saliva or secretions from the upper respiratory tract)
- or by air (airborne transmission by fine particles of less than 5 microns)

Comply with



OPPBTP specifications


MANDATORY in case the installer is working within one meter of another person:

- Wearing a type 1 non-sanitary mask (filtration greater than or equal to 90%)
- Wearing individual masks to be used by professionals in contact with the public, as FFP1 or as surgical type.



- Surgical mask, single use (duration 4 hours)
 - Great ease of breathing
 - Perfectly match the contours of the face thanks to its anatomical folding (3 folds)
 - High filtration efficiency
- EFB > 95% Bacterial Filtration Efficiency (mini)

- White internal face, Blue outer face.
- Hypoallergenic 3 layers, fiberglass free, latex free.
- Attachment by elastic loops, wide and comfortable
- Unique size
- Flat format 175 x 95mm



Packed in boxes of 50pcs



50 pcs



=



* Non contractual visuals